


## Easy-Shake Dispenser

Easy-Shake is a revolutionary hand held dispenser. It allows to dispense **one** applicator at a time to the **last** without touching other applicators.



The dispenser has a semi-cylindrical shape with a flattened area and a rounded area. The flat side facilitates guiding the applicators to a passageway. The shield covering separates the main body from the passageway. The purpose of separating the main body from the passageway is to prevent applicators to pile on top of each other in the passageway. As a result only one applicator is dispensed at a time – to the last applicator. Only one applicator is touched at a time as only one applicator is dispensed.

### Holding Position of Easy-Shake




Den-Brush is a Trademark of Denbur, Inc.  
Patented   
Manufactured by Denbur, Inc. ©2-08  
Made in USA

### Denbur, Inc.

P.O. Box 3473, Oak Brook, Illinois 60522-3473 USA  
630-986-9667 • 800-992-1399 • fax: 630-986-9688

E-mail: [denbur@denbur.com](mailto:denbur@denbur.com)  
[www.denbur.com](http://www.denbur.com)

# Den-Brush™ for Interdental Cleaning



## Disposable Triangular Fibrous Interdental Brush

Denbur, Inc.





**Den-Brush**

**Disposable Triangular Fibrous  
Interdental Brush**



The conventional interdental brushes sustain a thin wire with bristles cylindrically integrated within the wire. Such cleaning or flossing devices may be suitable for larger interdental areas. However, due to their cylindrical design they may irritate the gum line. Specifically because the bristles extend radically from the wire, the bristles form a generally cylindrical brush that does not correspond to the natural anatomy of the interdental spaces. Thus, the bristles frictionally interfere with and/or undesirably contact the gum line to irritate the sensitive gum areas.

**Den-Brush** is a new concept in interdental cleaning device. It is coated with bristles on 3 sides; the edge, and portion of both opposing side walls, while the back wall is void of any bristles. The back wall corresponds to or is positioned with respect to the user's gum line. A small hollowed area in the back wall prevents Den-Brush to place unnecessary pressure on the gum line. The back wall is void of bristles. As a result, there is no bristle contact between the back wall and the gum line. Den-Brush will not irritate the user's gum line.

There are two indentations in the handle portion to allow Den-Brush to bent into desired position. Both ends of the handle are rounded while the middle area consist of two flat sides for controlled and convenient holding position. Handle length allows the user to easily reach all areas of upper and lower jaw.



Coated tip fibrous bristles on 3 sides. Rounded end handle with two flat middle sides



Bent position for desired angle.

**Den-Brush/Den-Pick**

**938 939**  
Introductory Combo Pack  
6 Den-Brush/6 Den-Pick  
per dozen

**938 012**  
Introductory Pack  
12 Den-Brush per dozen

**939 012**  
Introductory Pack  
12 Den-Pick per dozen

**938 120**  
Economy Easy-Shake  
120 Den-Brush

**939 120**  
Economy Easy-Shake  
120 Den-Pick



Need photo  
with green



**Den-Brush** back wall or gum line void of fibers.



**Den-Brush** fibrous portion of both opposing side walls.



**Den-Brush** the fibrous edge.



**Den-Pick** is identical to Den-Brush structural body without the fibers.